



Certificate of Specification

GAL-G7 (KG7-100S)

| | |
|------------------------------------|---|
| CAS: | 96597-16-9 |
| Form: | White to yellowish powder |
| Formula: | $C_{50}H_{77}NO_{38}$ |
| Formula weight: | 1300.1 |
| Solubility (7%, H ₂ O): | Clear solution |
| EPS (HPLC.): | ≥ 90 % |
| H ₂ O (KF): | ≤ 3 % |
| 4-NP-G7 (enzym.): | ≤ 0.1 % |
| 4-Nitrophenol free: | ≤ 0.01 % |
| Storage: | 2 - 8 °C. Protect from light and moisture |
| Stability: | Stable at least 2 years at 2 - 8 °C. |