



KATALİZMED
İç ve Dış Ticaret San. Ltd. Şti.

Certificate of Specification

Cholesterol esterase (KCHE-101E)
Steryl Ester Acylhydrolase (EC 3.1.1.13)

Batch#

Expiry:

Specifications

Results

Origin :	Microorganism	Conform
Appearance :	White to light brown lyophilized powder	Conform
Activity :	≥ 25 u/mg	
Storage :	-20 °C	

The batches meet the range of specifications according to the above information.

Certifications are confirmed and signed according to the batch specifications.

Date: