



Certificate of Specification

Albumin Bovine (KALB-100)

Batch#

Expiry:

Specifications

CAS	9048-46-8
Molecular weight:	Approx. 66 kDa
Appearance	Off white to greyish powder or flakes
Purity:	> 98 %
pH:	6.5 - 7.5
Solubility in water (2%):	Total
Fat:	< 1 %
Total Nitrogen:	> 14 %
Total aerobic bacteria:	< 1000/g
E. coli:	Negative in 10 g
Salmonella:	Negative in 0.5g
Loss on drying:	< 6 %
Sulfated ashes	< 3 %
Storage:	2 - 8 °C recommended
Stability:	At least 2 years at 2 - 8 °C

The batches meet the range of specifications according to the above information.

Certifications are confirmed and signed according to the batch specifications.

Date: