



Certificate of Specification

HDAOS (KHDS-100)

Batch#

Expiry:

Specifications

CAS: 82692-88-4

Appearance : White powder

Molecular formula : $C_{11}H_{16}NO_6SNa$

Molecular weight: 313.30

Assay: > 99 %

Water content: < 0.5 %

Solubility (5%, H₂O): Clear, colorless solution

Storage: 2 – 8 °C. Protect from moisture and light.

Stability: At least 2 years at 2 - 8°C.

The batches meet the range of specifications according to the above information.

Certifications are confirmed and signed according to the batch specifications.

Date: