



Certificate of Specification

Arsenazo (KASN-100)

Batch#

Expiry:

Formula: $C_{22}H_{18}As_2O_{14}N_4S_2$

Molecular Weight: 776.37

Specifications

Appearance : Deep brown powder

Assay (HPLC) : $\geq 85 \%$

Specific extinction : 300 - 450

Sulfated ash : $\leq 0.2 \%$

Maximum of absorption : 535 - 543

Loss on drying : $\leq 8 \%$

Calcium : ≤ 100 ppm

Iron : ≤ 100 ppm

Copper : ≤ 50 ppm

Storage : Room Temperature

The batches meet the range of specifications according to the above information.

Certifications are confirmed and signed according to the batch specifications.

Date: