



KATALİZMED
İç ve Dış Ticaret San. Ltd. Şti.

Certificate of Specification

Ascorbate oxidase (KASCO-100E)
L-Ascorbate: oxygen oxidoreductase (EC 1.10.3.3)

Batch#

Expiry:

Specifications

Origin:	Cucurbita sp.
Appearance:	Light blue lyophilized powder
Activity:	≥ 200 u/mg
Contaminants:	Catalase: ≤ 0.1 % Phosphatase: ≤ 0.02 %
Storage:	-20 °C

The batches meet the range of specifications according to the above information.
Certifications are confirmed and signed according to the batch specifications.

Date: