



Certificate of Analysis

CKMM antibody (KCKM-100A)

Specifications

Description:	The KCKM-100A is a mixture of monoclonal antibodies
Appearance:	Deep-frozen solution
Immunogen:	Human Creatinine Kinase MM
Buffer:	10 mM Phosphate buffer, 9 g/L NaCl, 1g/L preservative
Specificity:	Human CK-MM: $\geq 99 \%$ Human CK-BB: $\leq 1 \%$
pH (25 °C):	7.4
Concentration:	≥ 17 mg/ml
Storage:	Store at -20 °C Can be stored at 2-8 °C for short term
Stability:	Stable for at least 24 months at -20 °C