

Certificate of Analysis

Lipoproteine lipase “Chemically Stabilized” (KLPL-100E)
Triacylglycero-protein acylhydrolase (EC 3.1.1.34)

Batch#

Expiry:

Specifications

Origin :	Pseudomonas sp.		
Appearance :	Light brown amorphous powder lyophilized		
Activity :	≥ 30 u/mg		
Contaminants :	Phosphatase :	≤ 10 ⁻³ %	
	Catalase :	≤ 2x10 ⁻² %	
	NADH oxidase :	≤ 10 ⁻³ %	
	Cholesterol oxidase :	≤ 2x10 ⁻³ %	
Storage :	-20 °C		

The batches meet the range of specifications according to the above information.

Certifications are confirmed and signed according to the batch specifications.

Date: