



Certificate of Specification

CRP Antibody polyclonal microorganism (KCRP-100A)

Batch#

Expiry:

Specifications

Source:

Goat

Specificity:

Human C-reactive protein (CRP)

Presentation:

liquid, 0.05M TRIS at pH 7.5, 0.5M NaCl, 15mM NaN₃

Assay (Becker Titer):

≥ 8.5 mg/ml

Total Protein (Biuret):

≥ 75 mg/ml

Storage:

store at 2 - 8 °C

The goats are inspected and certified free of Rinderpest, foot and mouth disease, contagious bovine pleuropneumonia and negative for Blue tongue.

The batches meet the range of specifications according to the above information.

Certifications are confirmed and signed according to the batch specifications.

Date: